

Incubators

PERFECTLY COORDINATED. PERFECTLY CONTROLLED.



2





Stable. Safe. Sensitive.

Memmert incubators for microbiology. Energy efficient, precise, 100% AtmoSAFE.

Even slight temperature deviations in the working chamber of an incubator may cause a test to fail. For this reason, the heating and control system of Memmert incubators are perfectly adapted to each other. During heating up and cooling down as well as in running operation, all appliances precisely keep the desired parameters within the smallest tolerance limits. Not only at one measuring point, but in the entire working chamber. Each individual Memmert incubator complies with the strict requirements of DIN 12880:2007-05 and is equipped with a maximum of safety functions. Each individual Memmert incubator is 100% AtmoSAFE.





INCUBATOR I	PAGE 4 TO 5
INCODATONI	I AGL 4 IO 3

TECHNICAL DATA PAGE 6 TO 7

Microbiological tests, colony counts, virology, toxicology

CO, INCUBATOR ICOmed **PAGE** 8 TO 9

TECHNICAL DATA PAGE 10 TO 11

Medical device class IIa for cultivation of cells or tissue, in-vitro fertilisation, gene expression

COMPRESSOR-COOLED INCUBATOR ICP PAGE 12 TO 13

TECHNICAL DATA PAGE 14 TO 15

Microbiological tests, colony counts, virology, toxicology, cultivation above and below room temperature, alternate stability tests

PELTIER-COOLED INCUBATOR IPP PAGE 16 TO 17

PAGE 18 TO 19

Protein crystallography, microbiological tests, colony counts, virology, toxicology, cultivation above and below room

temperature, alternate stability tests

COOLED STORAGE INCUBATOR IPS PAGE 20 TO 21

TECHNICAL DATA PAGE 22 TO 23

Microbiological tests, cultivation above and below room temperature, stability tests

SOFTWARE FOR MODELS U, UF TS, UNpa, S, I, ICOmed, ICP, IPP, IPS, HPP, ICH PAGE 24

AtmoCONTROL

TECHNICAL DATA

OPTIONS AND ACCESSORIES PAGE 25 TO 26

Available for all products


FEATURES OF MODEL VARIANTS

SingleDISPLAY and TwinDISPLAY



Incubator IN/INm and IF/IFm with SingleDISPLAY Incubator INplus/INmplus and IFplus/IFmplus with TwinDISPLAY Natural convection or forced air circulation AtmoCONTROL software

Model sizes: 30 / 55 / 75 / 110 / 160 / 260 / 450 / 750 +30 °C to +80 °C







As little air circulation as possible in the incubator


Forced air convection may destroy the protective layer from moist air that is generated during incubation over the samples. This would lead to dehydration of the culture. In a Memmert incubator, the perfect combination of all-round surface heating and temperature control system ensures that incubation generally takes place without forced air circulation. Provided the chamber is fully loaded and forced air circulation is required, it can be precisely adjusted in 10 % steps from 0 to 100 %.

Sterilisation

The chamber of the incubators INplus/IFplus/INmplus and IFmplus, including all installations and sensors, can be sterilised at +160 °C in a 4-hour programme to guarantee optimum hygiene.

Fresh air is preheated

Temperature deviations caused by fresh air can influence sample characteristics or prolong drying. In Memmert incubators, the fresh air is therefore fed through a pre-heating chamber and seamlessly introduced into the working chamber.



Air supply from outside



Memmert incubator Im is a medical device:

The incubator Im is a Class I medical device in accordance with the EU directive 93/42/EEC. In accordance with the intended use incubators INmplus and IFmplus may be used for warming of rinsing solutions and infusions as well as of contrast agents. Incubators IFm (with extended overtemperature protection — option A6) and IFmplus are used for warming non-sterile cloths and blankets. Incubators INm (with option A6) are intended for heating fango, silicate and APS packs for physical therapy and keeping them warm.

INCUBATORS I

according to DIN 12880:2007-05, EN 61010-1 (IEC 61010-1), 61010-2-010

Standard units are safety-approved and bear the test marks: (EAC not valid for medical devices)









Standard equipment

Interior: Stainless steel, material 1.4301 (ASTM 304) with all-round

deep-drawn ribs to integrate the large-area heating with

ceramic-metal sheath

Stainless steel grids, electropolished Internals:

(size 30 and 55: 1, size 75 to 750: 2)

Textured stainless steel, rear zinc-plated steel, intuitively Housing:

operated SingleDISPLAY or TwinDISPLAY (TFT colour display) with touchscreen; inner glass door, outside fully insulated

stainless steel door (from size 450 two leaves)

Fresh air: Admixture of pre-heated fresh air by electronically

adjusted airflap


Connection: Mains cable with plug (German type)

Installation: 4 feet; size 450 and 750 mounted on lockable castors

Interfaces:

LAN D





Model sizes/Description				30	55	75	110	160	260	450	750
Stainless steel interior	Volume		approx. I	32	53	74	108	161	256	449	749
	Width	(A)	mm	400	400	400	560	560	640	1040	1040
	Height	(B)	mm	320	400	560	480	720	800	720	1200
	Depth (less 39 mm for fan)	(C)	mm	250	330	330	400	400	500	600	600
	Stainless steel grids, electropolished (standard equipment)		number		1			2	2		
	Max. number of grids/shelves		number	3	4	6	5	8	9	8	14
	Max. loading per grid/shelf		kg			2	0			3	0
	Max. loading of chamber		kg	60	80	120	175	210		300	
Textured stainless steel	Width	(D)	mm	585	585	585	745	745	824	1224	1224
exterior	Height (size 450, 750 with castors)	(E)	mm	704	784	944	864	1104	1183	1247	1726
	Depth (without door handle), door handle + 56 mm	(F)	mm	434	514	514	584	584	684	784	784
Further data	Electrical load at 230 V, 50/60 Hz		approx. W	1600	1000	1250	1400	1600	1700	1800	2000
	Electrical load at 115 V, 50/60 Hz	approx. W	800			900			1500	1800	
	Working-temperature range		°C	min. 5 (IN/INplus/INm/INmplus) 10 (IF/IFplus/IFm/IFmplus) above ambient temperature up to +80							
	Setting temperature range		°C				+20 to	0 +80			
	Setting accuracy		°C				0.	.1			
Packing data	Net weight		approx. kg	48	57	66	76	96	110	161	217
	Gross weight (packed in carton)		approx. kg	64	76	85	101	122	161	227	288
	Width		approx. cm	66	73	73	83	83	93	133	133
	Height		approx. cm	89	95	113	105	130	138	144	191
	Depth		approx. cm	65	67	67	80	80	93	105	105
Order No. Incubators = Incubator				IN30 IN30m	IN55 IN55m	IN75 IN75m	IN110 IN110m	IN160 IN160m	IN260 IN260m	IN450 IN450m	IN750 IN750m
N = Natural convection F = Forced convection m = Medical device				IN30plus IN30mplus	IN55plus IN55mplus		IN110plus IN110mplus				
plus = Model with TwinDI	SPLAY			IF30 IF30m	IF55 IF55m	IF75 IF75m	IF110 IF110m	IF160 IF160m	IF260 IF260m	IF450 IF450m	IF750 IF750m
				IF30plus IF30mplus	IF55plus IF55mplus		IF110plus IF110mplus				

Options	30	55	75	110	160	260	450	750
Voltage 115 V, 50/60 Hz)	⟨2			
Extended overtemperature protection by additionally integrated Pt100 sensor for independent temperature monitoring for models IN/IF/INm/IFm				A	46			
Chamber modification for the appplication of reinforced perforated stainless steel shelves or stainless steel grids (bearing rails mounted in the working chamber) – includes replacement of 2 standard grids by 2 reinforced grids				_			K	1
Interior lighting for observing the load				F	R0			
Interior socket can only be ordered with limited temperature range up to max. +70 °C, ampacity 230 V/2.2 A, can be switched off with the On/Off switch, cannot be switched individually, moisture tight IP68 (only with SingleDISPLAY)				F	R3			
Entry port, 23 mm clear diameter, for introducing connections at the side, can be closed by flap and silicone stopper, standard positions left centre/centre left centre top right centre/centre right centre top					-0 -1 -2 -3			
Entry port, 23 mm clear diameter, can be closed by flap, in special positions (please, state location) left right rear				/////i	-4 -5 -6			
Entry port, 14 mm clear diameter, can be closed by flap, in special positions in the back wall (please, state location)				[06			
Entry port, 38 mm clear diameter, can be closed by flap, in special positions in the back wall (please, state location)				l	-7			
Entry port, 57 mm clear diameter, can be closed by flap, in special positions in the back wall (please, state location)					8			
Entry port, 100 mm clear diameter, can be closed by flap, in special positions in the back wall (please, state location)		-			F	:9		
4 – 20 mA current loop interface (0 to 90 °C ≜ 4 to 20 mA) Temperature controller, actual value Temperature of a Pt100 sensor positioned flexibly in chamber for external temperature monitoring (max. 1 SingleDISPLAY, max. 3 TwinDISPLAY)					/3 /6			
Fan speed monitoring with switching off the heating and with alarm in case of failure optional only for IFplus/IFmplus				\	/4			
Works calibration certificate for 3 temperatures: +37 °C, +52 °C, +70 °C Standard works calibration certificate (measuring point chamber centre) at +37 °C				D00	0126			

Accessories	30	55	75	110	160	260	450	750	
Stainless steel grid, electropolished (standard equipment)	E28884	E20	164	E20	165	E28891	E20	182	
Additional reinforced stainless steel grid, electropolished, max. loading 60 kg; from size 450 with guide bars and fixing screws (only in connection with option K1). Please consider max. loading of chamber		-		E29	767	E29766	B32	2190	
Perforated stainless steel shelf	B29727	B03	916	B00	325	B29725	B00	328	
Additional reinforced stainless steel shelf, max. loading 60 kg; with guide bars and fixing screws (only in connection with option K1). Please consider max. loading of chamber				_			B32	B32191	
Stainless steel slide-in drip tray, 15 mm rim (may affect the temperature distribution) — cannot be used in connection with option K1	E02070	E02	072	E02	.073	E29726	E02	E02075	
Max. loading per slide-in drip tray (kg)		1,5		3 4			8	8	
Stainless steel slide-in drip tray, 15 mm rim, with guide bars and fixing screws (can be used only in connection with option K1)				-		B32	2763		
Stainless steel bottom drip tray, 15 mm rim (may affect the temperature distribution) – cannot be used in connection with option K1	B04356	B04	B04358		1359	B29722	B04	1362	
Max. loading per bottom drip tray (kg)		1,5			3	4	8	8	
Stainless steel bottom drip tray, 15 mm rim (can be used only in connection with option K1)				-			B34	1055	
Wall bracket for wall mounting	B29755	B29756	B29757	B29758	B29759		///-///		
Guarantee extension by 1 year			GA1Q5				GA2Q5		




CO₂ Incubator ICOmed with TwinDISPLAY Software AtmoCONTROL

Model sizes: 50 / 105 / 150 / 240 + 18 °C to +50 °C Humidity 40 to 97 % rh CO_2 concentration 0 to 20 % O_2 concentration 1 to 20 %

CO₂ INCUBATOR ICOmed Safety at all times. When it comes to safety and user friendliness, the highly modern CO₂ incubator ICOmed is the perfect solution: Thanks to the battery-buffered ControlCOCKPIT, the operating display, logging and CO₂ control remain fully functional even when there is a power failure. All parameters are logged in accordance with the FDA and, when individually adjusted ranges for CO₂, O₂, temperature and humidity are exceeded, notifications can be sent to a mobile phone in addition to an alarm.

The control technology is so finely tuned that the setpoint temperature is guaranteed to be reached without temperature overshoots. With its rounded corners, the interior is easy to clean and can be sterilised for 60 minutes at 180 °C (including all sensors).

All ICOmed models are classified as medical product class IIa.







Comfort options for every application

These are only a few of the numerous features that don't let anything open to be desired:

- Two gas connections with quick release connectors for automatic switch-over of gas cylinders
- Electropolished working chamber
- Electronic control for active humidification and dehumidification (40 to 97 % rh)
- Control of oxygen concentration by introducing nitrogen, adjustment range from 1 % to 20 % O₃

Unrivalled user friendliness

All parameters can be set easily and intuitively both with the ControlCOCKPIT or the AtmoCONTROL software. The shutter box can be opened, allowing fast access to controls. Maintenance is possible even if the appliances are stacked. The appliance has USB and Ethernet connections as well as a data logger with a ten-year storage capacity. Data can be read and programmes can be transferred by remote access.

Minimising vaporisation and condensation

The active humidity control minimises vaporisation in the interior and ensures short recovery times after the door has been opened. Together with the heating of the interior from all six sides including the heated inner glass door, it prevents the dangerous formation of condensation and offers maximum protection for cell and tissue cultures. The turbulence-free chamber ventilation ensures a constant and uniform atmosphere.



CO, Incubator ICOmed is a medical device:

Memmert subjected its CO₂ incubator ICOmed to a comprehensive medical device evaluation. Every Memmert CO₂ incubator ICOmed is classified as a Class IIa medical device. The ICOmed is intended for the creation and maintenance of constant environmental conditions for application in the field of in vitro fertilisation (IVF), especially for the incubation of oocytes, spermatozoa and zygotes in special culture dishes for IVF application as well as for gene expression and the biosynthesis of RNA and proteins. The CE label on the appliances includes the mark 0197, denoting TRLP – TÜV Rheinland as the notified body.



CO, INCUBATORS ICOmed

with standard sterilisation programme (Humidity and ${\rm CO_2}$ sensor sterilised inside the ${\rm CO_2}$ incubator)

according to 12880:2007-05, EN 61010-1 (IEC 61010-1), 61010-2-010

Standard units are safety-approved and bear the test marks:







Standard equipment

Interior: Stainless steel, material 1.4301 (ASTM 304),

deep-drawn, seamlessly welded

Perforated stainless steel shelves Internals:

size 50: 1, sizes 105 - 240: 2;

and 1 stainless steel water dish (all sizes)

Textured stainless steel, rear zinc-plated steel, Housing:

intuitively operated TwinDISPLAY (TFT colour display) with touchsceen; fully insulated stainless steel door and heated inner glass door


Mains cable with plug (German type) Connection:

Installation: 4 adjustable feet

Interfaces:







Model sizes/Description	1		50	105	150	240
Stainless steel interior	Volume	approx. l	56	107	156	241
	Width (A)	mm	400	560	560	600
	Height (B)	mm	425	480	700	810
	Depth (less 35 mm for fan) (C)	mm	330	400	400	500
	Stainless steel shelves, perforated (standard equipment)	number	1	2	2	2
	Max. number of perforated shelves	number	5	6	10	12
	Max. loading per perforated shelf	kg		1	5	
	Max. loading of chamber	kg	75	90	120	140
tainless steel exterior	Width (D)	mm	559	719	719	759
	Height (variable through adjustable feet) (E)	mm	791	846	1066	1176
	Depth (without door handle, depth of door handle 56 mm) (F)	mm	521	591	591	691
	Fully insulated stainless steel door					
	Additional heated inner glass door					
urther data	Electrical load at 230/115 V, 50/60 Hz	approx. W	1000	1500	2000	2000
	Working-temperature range Standard sterilisation programme: 60 minutes at 180 °C (without removing the sensors)	°C	5 8	above ambient ter	mperature up to +	50
	Setting temperature range	°C		+18 t	o +50	
	Setting accuracy	°C		0	.1	
	Temperature fluctuations with time at 37 °C (to DIN 12880:2007-05)	K		+/-	0.1	
	Temperature variation in chamber at +37 °C (to DIN 12880:2007-05)	K		+/-	0.3	
	Humidity limitation thanks to a Peltier element; when water dish is full and inserted, the Peltier element limits the value of relative humidity in the interior to 93 % $$ rh +/- 2.5 %				+/- 0.3	
	Setting range active humidity control (with option K7)	% rh		40 to 97	and rh-Off	
	Setting accuracy humidity	% rh		0	.5	
	Digital electronic CO ₂ control with dual beam NDIR system, with auto-diagnostic system and acoustic fault indication, barometric pressure compensation	ty % rh ontrol with dual beam NDIR system, stem and acoustic fault indication,				
	Adjustment range CO ₂	% CO ₂		0 to	20	
	Variation in time CO ₂	% CO ₂		+/-	0.2	
	Setting accuracy CO ₂	% CO ₂		0	.1	
	Adjustment range O ₂	% O ₂		1 to	20	
	Setting accuracy O ₂	% O ₂		0	.1	
tandard accessories	Stainless steel water dish		1	1	1////1	1
	Works calibration certificate (measuring point chamber centre) at +37 °C, 5 % $\rm CO_2$ for standard units			[-	
	Works calibration certificate at 37 °C, 5 % $\rm CO_2$ and 90 % rh (requires option K7); standard equipment for units with active humidity control			1	-	
	Works calibration certificate at 37 °C, 5 % $\rm CO_2$, 90 % rh and 10 % $\rm O_2$ (requires option K7 and option T6); standard equipment for units with $\rm O_2$ control			1	-	

Model sizes/Descri	ption		50	105	150	240	
Packing data	Net weight	approx. kg	55	75	90	110	
J	Gross weight (packed in carton)	approx. kg	74	100	116	145	
	Width	Width	approx. cm	73	80	80	84
	Height	approx. cm	95	103	125	136	
	Depth	approx. cm	64	80	80	90	
Order No. CO ₂ Inc	cubators		ICO50med	ICO105med	ICO150med	ICO240med	

Options	50	105	150	240	
Voltage 115 V, 50/60 Hz		>	(2		
Battery-buffered ControlCOCKPIT Uninterrupted supply for the entire display unit (ControlCOCKPIT) and therefore complete documentation of all parameters even when there is a power failure. The CO, parameter is continuously regulated		(
Two gas connections with quick release connectors for automatic switch-over of gas cylinders					
Electropolished interior		///////////////////////////////////////	-2		
Active microprocessor control for humidifying and dehumidifying $(40 - 97 \% \text{ rh})$, incl. digital indication and auto-diagnostic system ensures even more rapid reaching of set humidity and very short recovery times while avoiding condensate formation. Humidity supply with distilled water (from an external tank) by a self-priming pump; integral bacteria block by generating hotsteam, dehumidifying via sterile filter		k	(7		
Control of oxygen concentration by N_2 inlet; adjustment range 1 % up to 20 % O_2 ; setting accuracy 0.1 %. (requires option K7)		T6			
Peltier cooling unit enables a working temperature of 37 °C even at higher ambient temperatures of up to 35 °C	K5				
Capacitive humidity sensor for measuring and displaying the relative humidity		k	(6		
Entry port (silicone), 40 mm clear diameter, for introducing connections, moisture tight, can be closed by silicone stopper, at the back, centre right; not available for ICO50med with active humidity control (option K7) or humidity display (option K6)		F	:7		
Inner door with partitioned glass doors	-		K4		
4 – 20 mA current loop interface $ \begin{array}{c} \text{Temperature controller, actual value (0 to +70 °C $\triangleq 4$ to 20 mA)} \\ \text{Humidity controller, actual value (0 – 100 % rh $\triangleq 4$ – 20 mA)} \\ \text{CO}_2 \text{ controller, actual value (0 – 25 % CO}_2 $\triangleq 4$ – 20 mA)} \\ \text{O}_2 \text{ controller, actual value (0 – 25 % O}_2 $\triangleq 4$ – 20 mA)} \\ \end{array} $		\ \	/3 /7 /9 /1		
Works calibration certificate for 5 %, 7 % and 10 % CO ₂ (measured at +37 °C) special works calibration certificates upon request		D00)106		
Start-up of ICOmed incubators and brief training (D, A, CH only), through Memmert service		k	(9		

Accessories	50	105	150	240
Additional perforated stainless steel shelf	E35160	E3	7418	E35158
Additional water dish	B38737		B38000	
Subframe (622 mm high) adjustable in height (sizes 150/240: should not be used for 2 stacked units)	B33504	В3	3505	B33506
Subframe (130 mm high); sizes 150/240: should not be used for 2 stacked units	B33507	В3	3508	B33509
HEPA-filter for chamber (filter class E11) according to EN 1822, packed in sterile condition, incl. fixing unit		В3	8739	
CO ₂ pressure reducing valve to DIN 8546, incl. gas cylinder monitor		E0	2087	
N, pressure reducing valve to DIN EN ISO 2503, incl. gas cylinder monitor (requires option T6)		E0	6162	
CO ₂ connection set, hose with coupling and clamp		В0	3881	
Central water supply, with filter cartridges for connection to the domestic water supply, only in combination with option K7. Product information on demand		ZV	/VR6	
Central water supply, without filter cartridges for connection to the domestic water supply (only for demineralised water in accordance with VDE 0510/DIN EN 50272), only in combination with option K7. Product information on demand		ZV	/VR7	
Guarantee extension by 1 year		G <i>A</i>	\3Q5	
Celltron benchtop shaker (accessories upon request)	-		E06724	



Compressor-cooled incubator ICP with TwinDISPLAY
AtmoCONTROL software

Model size: 55 0 °C to +60 °C

Model sizes: 110 / 260 / 450 / 750

-12 °C to +60 °C


COMPRESSOR-COOLED INCUBATOR ICP Ideal at temperatures around zero and below! If rapid and precise alternation between heating up and cooling down times in ramp operation is required, cooled incubators with compressor cooling prove to be in peak form — yet still work extraordinarily quiet. Due to the finely adjusted control technology, temperatures exactly reach the set point values without energy-intensive bursts of power.



Completely enclosed working chamber

Cooling and heating units are situated outside the working chamber inside the air jacket temperature control system surrounding the entire chamber interior ensuring quick and precise temperature control. The motor-driven forced air circulation, adjustable in 10 % steps via the ControlCOCKPIT ensures optimum temperature distribution.






ICP air jacket temperature control system

Integrated energy saving function

The cooling unit works extremely energy-efficiently, as there is no continuous heating against cooling. An intelligent DEFROST function enables defrosting as required.



COMPRESSOR-COOLED INCUBATORS ICP

according to DIN 12880:2007-05, EN 61010-1 (IEC 61010-1), 61010-2-010

Standard units are safety-approved and bear the test marks: $\mathbf{C} \in \mathbf{IH}$



Standard equipment

Stainless steel, material 1.4301 (ASTM 304) Interior:

Stainless steel grids, electropolished Internals:

(size 55: 1 grid, sizes 110 to 750: 2 grids)

Textured stainless steel, rear zinc-plated steel, intuitively Housing:

operated TwinDISPLAY (TFT colour display) with touchscreen;

fully insulated stainless steel door (from size 450 two leaves)


Connection: Mains cable with plug (German type)


Installation: mounted on lockable castors

Interfaces:









Max. number of grids Max. loading per grid Max. loading of cham Textured stainless steel exterior Width Height (with castors) Depth (without door left) Further data Electrical load at 230, Working-temperature (not suitable for long-	electropolished (standard equipment) s/shelves l/shelf nber handle), door handle + 56 mm	(A) (B) (C) (D) (E) (F)	approx. I mm mm mm number number kg kg mm mm	53 400 400 330 1 4 80 585 1153 514	108 560 480 400 5 20 150 745 1233	9 824 1552	449 1040 720 600 2 8 3 200 1224 1613	749 1040 1200 600 14 0
Height Depth (less 33 mm for Stainless steel grids, et Max. number of grids Max. loading per grid Max. loading per grid Max. loading of chames the Max of	electropolished (standard equipment) s/shelves l/shelf nber handle), door handle + 56 mm l/115 V, 50/60 Hz	(B) (C) (D) (E)	mm number number kg kg mm mm	400 330 1 4 80 585 1153	480 400 5 20 150 745 1233	800 500 9 824 1552	720 600 2 8 3 200 1224	1200 600 14 0
Depth (less 33 mm for Stainless steel grids, et Max. number of grids Max. loading per grid Max. loading of chame Textured stainless steel exterior Width Height (with castors) Depth (without door left) Depth (without door left) Depth (without door left) Setting temperature references to Setting temperature references as well as the set of t	electropolished (standard equipment) s/shelves l/shelf nber handle), door handle + 56 mm l/115 V, 50/60 Hz	(C) (D) (E)	mm number number kg kg mm mm	330 1 4 80 585 1153	400 5 20 150 745 1233	500 9 824 1552	8 3 200 1224	600 14 0
Stainless steel grids, e Max. number of grids Max. loading per grid Max. loading of cham Width Height (with castors) Depth (without door left) Further data Electrical load at 2300 Working-temperature (not suitable for long-During permanent op Setting temperature references	electropolished (standard equipment) s/shelves l/shelf nber handle), door handle + 56 mm l/115 V, 50/60 Hz	(D) (E)	number number kg kg mm mm	1 4 80 585 1153	5 20 150 745 1233	9 824 1552	8 3 200 1224	14 0 1224
Max. number of grids Max. loading per grid Max. loading of cham Width Height (with castors) Depth (without door left) Further data Electrical load at 230, Working-temperature (not suitable for long- During permanent op Setting temperature ref	handle), door handle + 56 mm //115 V, 50/60 Hz	(E)	number kg kg mm mm	80 585 1153	20 150 745 1233	9 824 1552	8 200 1224	1224
Max. loading per grid Max. loading of cham Width Height (with castors) Depth (without door l Further data Electrical load at 230 Working-temperature (not suitable for long- During permanent op Setting temperature re	handle), door handle + 56 mm /115 V, 50/60 Hz	(E)	kg kg mm mm	80 585 1153	20 150 745 1233	824 1552	200 1224	1224
Max. loading of cham Textured stainless steel exterior Width Height (with castors) Depth (without door land) Further data Electrical load at 230, Working-temperature (not suitable for long- During permanent op Setting temperature re-	handle), door handle + 56 mm /115 V, 50/60 Hz	(E)	kg mm mm	585 1153	150 745 1233	1552	200 1224	1224
Textured stainless steel exterior Height (with castors) Depth (without door labeled at 230.) Working-temperature (not suitable for long-During permanent op Setting temperature recommended)	handle), door handle + 56 mm /115 V, 50/60 Hz	(E)	mm mm	585 1153	745 1233	1552	1224	
Further data Height (with castors) Depth (without door later than the properties of the properties o	/115 V, 50/60 Hz	(E)	mm	1153	1233	1552		
Height (with castors) Depth (without door left) Further data Electrical load at 230, Working-temperature (not suitable for long- During permanent op Setting temperature ref	/115 V, 50/60 Hz						1613	1950
Further data Electrical load at 230, Working-temperature (not suitable for long- During permanent op Setting temperature re	/115 V, 50/60 Hz	(F)	mm	514	F0.4			
Working-temperature (not suitable for long- During permanent op Setting temperature r	e range				584	33 1552 1613	784	
(not suitable for long- During permanent op Setting temperature r	e range		approx. W			1200		
	peration, the glass door may ice over)		°C		-12 to +	60 (ICP 55 () to +60)	
Setting accuracy	range		°C		-12 to +6	50 (ICP 55 -	5 to +60)	
			°C			0.1		
Packing data Net weight			approx. kg	89	113	157	217	249
Gross weight (packed	l in carton)		approx. kg	111	141	214	282	319
Width			approx. cm	76	88	93	133	133
Height			approx. cm	133	141	+60 (ICP 55 -5 to +60) 0.1 157 214 282	170	215
Depth			approx. cm	68	81	93	105	105
Order No. Compressor-Cooled Incubators				ICP55	ICP110	ICP260	ICP450	ICP75

Options	55	110	260	450	750
Voltage 115 V, 50/60 Hz			X2		
Chamber modification for the application of reinforced perforated stainless steel shelves or stainless steel grids (bearing rails mounted in the working chamber) – includes replacement of 2 standard grids by 2 reinforced grids		-		K	(1
Interior socket, ampacity 230 V/2.2 A, can be switched off with the On/Off switch, cannot be switched individually, moisture tight IP68			R3		
Entry port, 23 mm clear diameter, for introducing connections at the side, can be closed by flap, standard positions left centre/centre left centre top right centre top		_	F0 F1	F3	
Entry port (silicone), 40 mm clear diameter, for introducing connections, moisture tight, can be closed by silicone stopper, at the back (please, state location)			F7		
4 – 20 mA current loop interface (-20 to +70 °C			V3 V6		
Fan speed monitoring with switching off the heating and with alarm in case of failure			V4		
Works calibration certificate for 3 temperatures: 0 °C, +37 °C, +60 °C Standard works calibration certificate (measuring point chamber centre) at +10 °C and +37 °C			D00130		

Accessories	55	110	260	450	750
Stainless steel grid, electropolished (standard equipment)	E20164	E20165	E28891	E20	182
Additional reinforced stainless steel grid, electropolished, max. loading 60 kg; from size 450 with guide bars and fixing screws (only in connection with option K1). Please consider max. loading of chamber	_	E29767	E29766	B32	190
Perforated stainless steel shelf	B03916	B00325	B29725	B00	328
Additional reinforced stainless steel shelf, max. loading 60 kg; with guide bars and fixing screws (only in connection with option K1). Please consider max. loading of chamber		-		B32	191
Stainless steel slide-in drip tray, 15 mm rim (may affect the temperature distribution) — cannot be used in connection with option K1	E02072	E02073	E29726	E02	075
Max. loading per slide-in drip tray (kg)	1.5	3	4	8	3
Stainless steel slide-in drip tray, 15 mm rim, with guide bars and fixing screws (can be used only in connection with option K1)		/// - ///		B32	763
Stainless steel bottom drip tray, 15 mm rim (may affect the temperature distribution) — cannot be used in connection with option K1	B04358	B04359	B29722	B04	362
Max. loading per bottom drip tray (kg)	1.5	3	4	8	3
Stainless steel bottom drip tray, 15 mm rim (can be used only in connection with option K1)		<u>-</u>		B34	055



Peltier-cooled incubator IPP with SingleDISPLAY Peltier-cooled incubator IPPplus with TwinDISPLAY AtmoCONTROL software

Model sizes: 30 / 55 / 110 / 260 / 750

 $0 \, ^{\circ}\text{C} \text{ to } +70 \, ^{\circ}\text{C}$

PELTIER-COOLED INCUBATOR IPP Heating and cooling seamlessly with one system thanks to Peltier technology. In this respect, cooled incubators IPP not only contribute to climate protection, but it also achieves an additional decrease in operating costs of up 90 % compared to compressor technology. This perfect development from the environmentally friendly and energy-saving heating/cooling technology by Memmert convinces by outstanding control precision and extremely small fluctuations.





Extremely quiet and vibration-free


The fact that no compressor is required saves space and brings peace and quiet to the laboratory. As Peltier-cooled incubators IPP are almost vibration-free, they can also be applied in entomology. If defined humidity is also required, an alternative would be the constant climate chamber HPP, which is also equipped with Peltier technology.

No condensation in the interior chamber

Due to the closed Peltier cooling system, no outside air is exchanged. Physically derived, unavoidable formation of condensation during the cooling process does not take place in the interior chamber but on the outside heat sink. In addition, the in the Peltier elements integrated fans ensure a rapid transport of energy as well as an optimal temperature distribution.

Energy-saving heating/cooling technology combination

In contrast to compressor systems, Peltier technology is particularly economical at temperatures close to the ambient temperature, since energy is only required during heating or cooling. Therefore heating and cooling function are particularly precisely adjusted to each other.





PELTIER-COOLED INCUBATORS IPP

according to DIN 12880:2007-05, EN 61010-1 (IEC 61010-1), 61010-2-010

Standard units are safety-approved and bear the test marks: CE III



Standard equipment

Stainless steel, material 1.4301 (ASTM 304), deep-drawn Interior:

Stainless steel grids, electropolished Internals:

(sizes 30 and 55: 1, sizes 110 to 750: 2)

Textured stainless steel, rear zinc-plated steel, intuitively Housing:

operated SingleDISPLAY or TwinDISPLAY (TFT colour display)

with touchscreen

Double doors: Outside stainless steel, fully insulated, inside glass


(size 750 two-leaves)

Connection: Mains cable with plug (German type) Installation: 4 feet; size 750 mounted on lockable castors

Interfaces:

Ethernet T LAN D





Number of Peltier elements in the rear: Sizes 30 to 55: 1 Size 110: 2 Size 260: 3 Size 750: 6

Order No. Peltier-Cool			IPP30	IPP55	IPP110	IPP260	IPP750			
	Depth		approx. cm	65	67	80		105		
	Height		approx. cm	89	95	105	824 1183 774 820 temperature) to 0 114 165 93 138 93 IPP260	191		
	Width		approx. cm	66	73	83	93	133		
J	Gross weight (packed in carton)		approx. kg	56	71	103		301		
Packing data	Net weight		approx. kg	40	52	78	temperature) to 40 0 0 114 165	230		
	Setting accuracy		°C			0.1				
	Setting temperature range		°C			0 to +70				
	Working-temperature range without light Working-temperature range with light		°C °C	0 (at least 20 bel		mperature) to +	-70		
Further data	Electrical load at 230/115 V, 50/60 Hz		approx. W	140	275	550	820	1100		
	Depth (without door handle), door handle + 56 mm	(F)	mm	524	604	674	9 640 9 800 0 500 2 9 0 20 5 824 4 1183 4 774 0 820 ent temperature) to + +40 -70 114 8 165 93 5 138 93 10 IPP260	874		
exterior	Height (size 750 with castors)	(E)	mm	704	784	480 8 400 5 5 20 150 745 8 864 11 674 7 550 8 low ambient temperature 10 to +40 0 to +70 0.1 78 1 103 1 83 9 105 1 80 9 IPP110 IPP	1183	1726		
Textured stainless steel	Width	(D)	mm	585	585	745	824	1224		
	Max. loading of chamber		kg	60	80	150	20	00		
	Max. loading per grid/shelf		kg		2	20		30		
	Max. number of grids/shelves		number	3	4	5	9	14		
	Stainless steel grids, electropolished (standard equipment)		number		1		2			
	Depth (less 10 mm for fan – Peltier)	(C)	mm	250	330		640 800 500 2 9 824 1183 774 820 emperature) to 114 165 93 138 93 IPP260	600		
	Height	(B)	mm	320	400			1200		
stanness steer interior	Width	(A)	approx. I	400	400			1040		
Stainless steel interior	Volume		anney I	32	53	108	256	749		

Options	30	55	110	260	750	
Voltage 115 V, 50/60 Hz			X2			
Chamber modification for the application of reinforced perforated stainless steel shelves or stainless steel grids (bearing rails mounted in the working chamber) – includes replacement of 2 standard grids by 2 reinforced grids			-		K1	
Light module cold white 6,500 K: light strips arranged on the side walls of the interior, 10 strips for model 110, 14 for model 260/750, programme-controlled dimming from 0 to 100 % (in 1 % steps), ramp programming in combination with temperature (only with TwinDISPLAY)		_		T7		
Light module cold white 6,500 K + warm white 2,700 K: LED light strips – 10 strips for model 110, 14 for models 260/750 – (5 resp. 7 alternating cold white light strips and 5 resp. 7 warm white light strips) on the side walls of the interior, programme-controlled dimming from 0 to 100 % (in 1 % steps), ramp programming in combination with temperature (only with TwinDISPLAY)		_		Т8		
Light module warm white 2,700 K: light strips arranged on the side walls of the interior, 10 strips for model 110, 14 for model 260/750, programme-controlled dimming from 0 to 100 % (in 1 % steps), ramp programming in combination with temperature (only with TwinDISPLAY)		_		Т9		
Interior socket, ampacity 230 V/2.2 A, can be switched off with the On/Off switch, cannot be switched individually, moisture tight IP68			R3			
Entry port, 23 mm clear diameter, for introducing connections, can be closed by flap, standard positions (F0 and F2 not for model size 260 with light module; F0 – F3 not for model size 110 with light module) right centre/centre right centre top			F0 F1 F2 F3			
Entry port, 23 mm clear diameter for introducing connections, can be closed by flap (please, state location) left right rear			F4 F5 F6			
Entry port, 14 mm clear diameter, can be closed by flap, in special positions in the back wall (please, state location)	D6					
Entry port, 38 mm clear diameter, can be closed by flap, in special positions in the back wall (please, state location)	F7					
4 – 20 mA current loop interface (-10 to +80 °C ≙ 4 – 20 mA) Temperature controller, actual value Temperature of a Pt100 sensor positioned flexibly in chamber for external temperature monitoring (max. 1 SingleDISPLAY, max. 3 TwinDISPLAY)			V3 V6			
Works calibration certificate for 3 temperatures: +5 °C, +37 °C, +60 °C Standard works calibration certificate (measuring point chamber centre) at +10 °C and +37 °C			D00129			

Accessories	30	55	110	260	750
Stainless steel grid, electropolished (standard equipment)	E28884	E20164	E20165	E28891	E2018
Additional reinforced stainless steel grid, electropolished, max. loading 60 kg; size 750 with guide bars and fixing screws (only in connection with option K1). Please consider max. loading of chamber		_	E29767	E29766	B3219
Perforated stainless steel shelf	B29727	B03916	B00325	B29725	B0032
Additional reinforced stainless steel shelf, max. loading 60 kg; with guide bars and fixing screws (only in connection with option K1). Please consider max. loading of chamber			-		B3219
Stainless steel slide-in drip tray, 15 mm rim (may affect the temperature distribution) – cannot be used in connection with option K1	E02070	E02072	E02073	E29726	E0207
Max. loading per slide-in drip tray (kg)	1.	.5	3	4	8
Stainless steel slide-in drip tray, 15 mm rim, with guide bars and fixing screws (can be used only in connection with option K1)		-	-		B3276
Stainless steel bottom drip tray, 15 mm rim (may affect the temperature distribution) – cannot be used in connection with option K1	B04356	B04358	B04359	B29722	B0436
Max. loading per bottom drip tray (kg)	1.	.5	3	4	8
Stainless steel bottom drip tray, 15 mm rim (can be used only in connection with option K1)			-		B3405
Guarantee extension by 1 year		GA1Q5		ĞA.	2Q5



Cooled storage incubator IPS with SingleDISPLAY AtmoCONTROL software

Model sizes: 260 / 750 +14 °C to +45 °C

COOLED STORAGE INCUBATOR IPS Save energy and reduce the strain on the climate at the same time! If microbiological cultures, BOB5 samples, drinks containers or cosmetics need to be stored over a long period at constant temperatures, cooled storage incubators IPS with energy-efficient Peltier technology are the perfect choice: absolute reliability, precision, durability and eco-friendliness.





Considerable potential for savings in acquisition and operating costs

Temperature changes are not always necessary for long-term storage or incubating. So why design heating, cooling and controlling systems for rapid heating up and cooling down times? The performance of the IPS was tailor-made for permanent operation at constant temperatures close to room temperature. The advantage: Acquisition costs and operating costs are considerably reduced in comparison to conventional cooled incubators with compressor technology, as well as to a large Peltier-cooled incubator.



Ideal for high ambient temperatures

Thanks to Peltier elements integrated for cooling the working chamber, the chamber load won't break into sweat even at high ambient temperatures. Constant and precise incubation at room temperature is guaranteed.



Low in vibration and durable for absolutely safe long-term storage

Like the cooled incubator IPP, the IPS offers all the advantages of Peltier technology to the user. Its interior chamber that is completely insulated from the environment minimises the risk of drying out of the samples. It is practically noise-free and not only reduces stress on the chamber load but also soothes the nerves of employees thanks to its quiet operation.



Glimpse into a Memmert storage incubator: Peltier elements guarantee perfect climate inside the chamber.

COOLED STORAGE INCUBATORS IPS

according to DIN 12880:2007-05, EN 61010-1 (IEC 61010-1), 61010-2-010

Standard units are safety-approved and bear the test marks: CE



Standard equipment

Interior: Stainless steel, material 1.4301 (ASTM 304),

deep-drawn

2 stainless steel grids, electropolished Internals:

Textured stainless steel, rear zinc-plated steel, Housing:

intuitively operated SingleDISPLAY (TFT colour display) with touchscreen

Outside stainless steel, fully insulated, inside glass Double doors:


(size 750 two leaves)


Mains cable with plug (German type) Connection:

Installation: 4 feet; size 750 mounted on lockable castors

Interfaces:

Ethernet LAN D





2 Peltier elements in the rear

Model sizes/Descriptio	n			260	750
Stainless steel interior	Volume		approx. I	256	749
	Width	(A)	mm	640	1040
	Height	(B)	mm	800	1200
	Depth (less 10 mm for fan – Peltier)	(C)	mm	500	600
	Stainless steel grids, electropolished (standard equipment)		number	7////	2
	Max. number of grids/shelves		number	9	14
	Max. loading per grid/shelf		kg	20	30
	Max. loading of chamber		kg	20	00
Textured stainless steel exterior	Width	(D)	mm	824	1224
	Height (size 750 with castors)	(E)	mm	1183	1726
	Depth (without door handle), door handle + 56 mm	(F)	mm	774	874
Further data	Electrical load at 230/115 V, 50/60 Hz		approx. W	55	50
	Working-temperature range/Setting temperature range		°C	+14 to	0 +45
	Setting accuracy		°C	0.	.1
Packing data	Net weight		approx. kg	113	230
	Gross weight (packed in carton)		approx. kg	164	301
	Width		approx. cm	93	133
	Height		approx. cm	138	191
	Depth		approx. cm	93	105
out and contribution	and the beautiful and the second seco				

Order No. Cooled storage Incubators IPS260 IPS750

Options	260	750
Voltage 115 V, 50/60 Hz	XZ	2
Chamber modification for the application of reinforced perforated stainless steel shelves or stainless steel grids (bearing rails mounted in the working chamber) — includes replacement of 2 standard grids by 2 reinforced grids	-	K1
Interior socket, ampacity 230 V/2.2A, can be switched off with the On/Off switch, cannot be switched individually, moisture tight IP68	R3	3
Entry port, 23 mm clear diameter, for introducing connections at the side, can be closed by flap, standard positions left centre top left centre top right centre/centre right centre top	F1 F2	
Entry port, 23 mm clear diameter, can be closed by flap, in special positions (please, state location) lef righ rea	t F5	j
Entry port, 14 mm clear diameter, can be closed by flap, in special positions in the back wall (please, state location)	D6	5
Entry port, 38 mm clear diameter, can be closed by flap, in special positions in the back wall (please, state location)	F7	
4 – 20 mA current loop interface (0 to +70 °C ≜ 4 – 20 mA) Temperature controller, actual value Temperature of a Pt100 sensor positioned flexibly in chamber for external temperature monitoring (SingleDISPLAY) Works calibration certificate for a freely selectable temperature value Standard works calibration certificate (measuring point chamber centre) at +18 °C and +25 °C		5

Accessories	260	750
Stainless steel grid, electropolished (standard equipment)	E28891	E20182
Additional reinforced stainless steel grid, electropolished, max. loading 60 kg; size 750 with guide bars and fixing screws (only in connection with option K1). Please consider max. loading of chamber	E29766	B32190
Perforated stainless steel shelf	B29725	B00328
Additional reinforced stainless steel shelf, max. loading 60 kg; with guide bars and fixing screws (only in connection with option K1). Please consider max. loading of chamber	////-	B32191
Stainless steel slide-in drip tray, 15 mm rim (may affect the temperature distribution) — cannot be used in connection with option K1	E29726	E02075
Max. loading per slide-in drip tray (kg)	4	8
Stainless steel slide-in drip tray, 15 mm rim, with guide bars and fixing screws (can be used only in connection with option K1)	-	B32763
Stainless steel bottom drip tray, 15 mm rim (may affect the temperature distribution) – cannot be used in connection with option K1	B29722	B04362
Max. loading per bottom drip tray (kg)	4	8
Stainless steel bottom drip tray, 15 mm rim (can be used only in connection with option K1)	/////-	B34055
Guarantee extension by 1 year	GA	2Q5

SOFTWARE AtmoCONTROL

AtmoCONTROL – The innovative control and logging software

Parameters such as temperature and humidity as well as the process time can be set directly at the ControlCOCKPIT.

Ramp programming is done via the control and logging software AtmoCONTROL, which features a completely new software design.

Drag, drop & go!

Numerical and graphic programming of complex processes is a thing of the past. Today, programming is done via AtmoCONTROL by means of the mouse or touchpad on your notebook. Even the most complex ramp programmes are created within minutes. Simply drag & drop the graphical symbols for the desired parameters to the input field and change the values according to your wishes with a mouse click.

Programme functions SingleDISPLAY and TwinDISPLAY

- Reading out, managing and organising the data logger
- Saving the log memory in various formats
- Online monitoring of up to 32 connected appliances
- Optical alarms when the alarm limits individually set at the ControlCOCKPIT are exceeded
- Automatic alarm to one or several e-mail addresses



Additional functions TwinDISPLAY

- Intuitive programming and archiving of ramps and programme sequences
- Synchronous visualisation of the created programme sequence during programming
- Application-specific repeat functions (loops) can be inserted within a temperature control programme in any place
- · Simple creation of repeating weekly programmes
- Programming, managing and transferring programmes via Ethernet interface or USB port

SPECIAL EQUIPMENT FOR MODELS U, UF TS, UNpa, S, I, I	CP, IPP, I	PS, HF	PP, ICH								ICOmed
Options for models U, UF TS, UNpa, S, I, ICP, IPP, IPS, HPP, ICH	30 55	5 75	110	160	260	400	450	750	1060	1400	50 / 105 / 150 / 24
Door with lock (safety lock); for models UFTS per side; standard with SN/SF and SNplus/SFplus 450 and 750 (not for models ICOmed)				В	6						////-///////
Door hinged on the left; for models UF TS per side			B8					//-/			B8
Potential-free contact (24 V/2 A) with socket to NAMUR NE 28 for external monitoring (indicates when setpoint is reached); models ICOmed: when set points of temperature and CO ₂ are reached						Н	5				
Potential-free contact for combination error message (e.g. supply failure, sensor fault, fuse)						Н	6				
Potential-free contact (24 V/2 A) with socket to NAMUR NE 28, for signal generation, controlled by programme segment, for free-selectable functions to be activated (e.g. activation of audible and visual signals, exhaust motors, fans, stirrers, etc.). Only for units with TwinDISPLAY; max. 2 contacts on 1-phase appliances; max. 4 contacts on 3-phase appliances (not for models ICOmed) 2 contacts 4 contacts		H72 - H74							_		
Process-dependent door lock (only for units with TwinDISPLAY); for models UF TS see page 11 of oven brochure; not for models ICOmed				D	4						-
Door-open-recognition, incl. alarm, shuts down fan and after 30 sec. also heating (only for units with TwinDISPLAY); for models UF TS per side; standard with ICOmed, ICH C, ICH L				٧	5						_
Flexible Pt100 for positioning in chamber or in load with socket, 4-pin, according to NAMUR NE 28, for external temperature recording (load temperature) max. 3 sensors; not for models ICOmed	H4						Н4				-
Flexible Pt100 temperature sensor, positioned flexibly in chamber or load, for local temperature measurement (up to 3 additional sensors are possible). The measured temperature can, if required, be indicated on the display, recorded in the integral data store, and can be documented via the AtmoCONTROL software. Not for models ICOmed	H8				H8					<u>-</u>	
MobileALERT, notification by SMS in case of any error or alarm of the device. Requires option H6 "floating contact for alarm"		C3									
MobileALERT for up to 4 alarm notifications; standard: temperature and CO ₂ alarm, additionally humidity alarm (when equipped with option K7) and O ₂ alarm (when equipped with option T6)					-						C4
Temperature restriction (for UN/UF/UNplus/UFplus/UNm/UFm/UNmplus/ UFmplus and models UF TS); Temperatures: +60, +70, +80, +95, +100, +120, +160, +180, +200, +220 or +250 °C (Please, indicate upon ordering)			A8			-		A8			_
Castor frame (2-part), height 140 mm (not for models UF TS, ICP, ICH, ICH L, ICH C, ICOmed)			R9								

SPECIAL EQUIPMENT FOR MODELS U, UF TS, UN	pa, S, I	, ICP,	IPP,	IPS, HP	P, IC	.Н						ICU	med
Accessories for models U, UF TS, UNpa, S, I, ICP, IPP, IPS, HPP, ICH	30	55	75	110	160	260	400	450	750	1060	1400	50 / 105	/ 150 / 24
JSB-Ethernet adapter							E06	5192					
Ethernet connection cable 5 m for computer interface							E06	5189					
JSB User-ID stick (with User-ID licence): Oven-linked authorisation icence (User-ID-programme) on Memory-stick, prevents undesired manipulation by unauthorised third parties. When reordering please							B3:	3170					
specify serial number (only for units with TwinDISPLAY) JSB stick with documentation software AtmoCONTROL and operation manual for products with SingleDISPLAY							B3:	3171					
the standard equipment of appliances with TwinDISPLAY ncludes one USB stick with AtmoCONTROL) Set of height adjustable feet (4 pcs) not available for ICP, ICH,													
CH L, ICH C – standard on models ICOmed Stacking set (4 pcs) for stacking of appliances of same size			B2	9768									
not for models 160, 260, 450, 750, 1060, 1400, ICH110, ICH110L, CH110C, ICP55, ICP110) Stacking set (consisting of stacking corners, one connecting plate		B297	744					-				B29744	<u> </u>
for the rear, two wall brackets) for stacking two units ICO150med or ICO240med Plug-in tube extension (outer diam. 60.3 mm, inner 57 mm),							-						B42114 (1 B42115 (2
straight, for exhaust air ducting (if necessary for connection by hose), only models U, I, S not for models UF TS					B2	29718						-	
Plug-in tube extension (outer diam. 60.3 mm, inner 57 mm), angled, for exhaust air ducting (if necessary for connection by hose), only models U, I, S not for models UF TS					B2	29719						-	
Flush-fit unit (stainless steel frame covering gap between oven and wall opening), with air slots	B29728	B29730	B2973	2 B29734	B2973	6 B29738	B42116	B29740	B29742	B4211	3	-	
Flush-fit unit (stainless steel frame covering gap between oven and wall opening), without air slots for models UFTS see page 11 of oven brochure; not for models ICOmed	B29729	B29731	B2973	3 B29735	B2973	7 B29739	B42117	B29741	B29743	B4211)	-	
Subframe, adjustable in height size 30 to 75: height 600 mm, size 110 to 450: height 500 mm); not for models ICOmed, UF TS and HPP400	B29745	B297	747	B297	749	B29751	-	B29753	3		-		
ubframe, on castors size 30 to 75: height 660 mm, size 110 to 160: height 560 mm); oot for models ICOmed and UFTS	B29746	B297	748	B297	750					-			
Subframe, adjustable in height, height 130 mm, or example for units with fresh air filter; not for models ICOmed and UF TS	B33657	B336	659	B336	561	B33664					-		
Software conforming to FDA AtmoCONTROL. Meets the requirements for the use of electronically stored data sets and electronic signatures as laid down in Regulation 21 CFR Part 11 of the US Food and Drug Administration (FDA). Base licence for the control of one unit only for units with TwinDISPLAY)	s s FDAQ1												
ntegration of additional units (up to max. 15 units) into an already existent FDA-software licence (only for units with TwinDISPLAY)							FD	AQ2					
Q document with device-specific works test data, DQ/PQ check list as support for validation by customer Q/OQ document with device-specific works test data for one free-							D0	0124					
selectable temperature value, incl. temperature distribution survey at Memmert for 9 measuring points (size 30), 27 measuring points (sizes 55 – 1060) to DIN 12880:2007-05, PQ check list as support for validation by customer. Further temperature values and validation at customer site on demand	D00125					D00127						-	
Q/OQ document with device-specific works test data for one free- selectable temperature and humidity value, incl. temperature distri- bution survey at Memmert for 27 measuring points (26 measuring points on model HPP1400) to DIN 12880:2007-05, PQ check list as support for validation by customer (models HPP and ICH). Validation at customer site on demand		-		D00136	-	D00	136	-	D00136	-	D00136		_
Q/OQ document with device-specific works test data for one ree-selectable temperature, humidity and light value, incl. tempeature distribution survey at Memmert for 27 measuring points o DIN 12880:2007-05, PQ check list as support for validation by customer (models HPP with light and ICH L). Validation at sustomer site on demand		_		D00137	_	D00	137	_	D00137			-	
Q/OQ document with device-specific works test data for one ree-selectable CO ₂ , humidity and temp. value, incl. temp. distribution urvey at Memmert for 27 measuring points to DIN 12880:2007-05, Q check list as support for validation by customer (models ICH C nd ICOmed, on models ICOmed a free-selectable humidity value is nly possible with option K7). Validation at customer site on demand		<u>-</u>		D38897	-	D38897			D38897			D3	8897
Q/OQ document with device-specific works test data for one ree-selectable CO ₂ and temperature value, incl. temperature listribution survey at Memmert for 27 measuring points to DIN 12880:2007-05, PQ check list as support for validation by ustomer (model ICOmed). Validation at customer site on demand	– D388						8898						
xternal measuring instrument with sensors for daylight and UV-light. Product information on demand (models HPP, ICH L, IPPplus)					B0471	3				-	B04713		-
oitto with additional measuring head for temperature and umidity measurement. Product information on demand models HPP, ICH L, IPPplus)					B0471	4				-	B04714		_

SPECIAL EQUIPMENT FOR MODELS VO, VOcool, HCP, TTC, CTC							
Options for models VO, VOcool, HCP, TTC, CTC	200	400	500	108	153	246	256
Interface Ethernet instead of USB including software				W4			
RS232 interface instead of USB				W6			
Computer interface RS485 (for networking a max. of 16 ovens) instead of RS232				V2			
Door with lock (safety lock, not available for VO, VOcool, TTC/CTC)				В6			
Flexible Pt100 for positioning in chamber or in load with socket, 4-pin, according to NAMUR NE 28, for external temperature recording (load temperature); for VO and VOcool on demand	Н4						
Additional Pt100 temperature sensor, positioned flexibly in chamber or load, for local temperature measurement (up to 3 additional sensors are possible). The measured temperature can, if required, be indicated on the multifunction display, recorded in the integral ring store, and can be documented via the "Celsius" software or on an attached printer. (Not available for VO, VOcool, TTC and CTC)	Н8						
Potential-free contact (24 V/2 A) with socket, according to NAMUR NE 28 for external monitoring (indicates when setpoint is reached)	Н5						
Potential-free contact (24 V/2 A), with socket, according to NAMUR NE 28 for combination error message (e.g. supply failure, sensor fault, fuse)	H6						
Potential-free contact (24 V/2 A), with socket, according to NAMUR NE 28, triple, for signal generation, controlled by programme segment for a total of 3 freely selected functions to be activated (e.g. acoustic and visual signals, exhaust motors, fans, stirrers etc.). Not available for VO, VOcool				Н7			
MobileALERT, notification by SMS in case of any error or alarm of the device. Requires option H6 "floating contact for alarm"				C3			

Accessories for models VO, VOcool, HCP, TTC, CTC	200	400	500	108	153	246	256				
USB connection cable for computer interface				E03643							
Parallel/USB converter cable with integrated power supply unit to connect HP printers with USB interface to MEMMERT units				E05300							
Documentation package consisting of parallel USB converter cable including PCL3-compatible HP colour inkjet printer with USB interface (HP OfficeJet 6000 or successor) for direct connection of printer to Memmert unit				B04432							
Temperature profile write/read unit for programming via PC, for writing to and reading from the chip card, up to 40 ramps				E05284							
Additional chip card, blank, formatted (32 kB MEMoryCard XL for a maximum of 40 ramps)				E04004							
Oven-linked authorisation card (User-ID-Card) prevents undesired manipulation by unauthorised third parties. When reordering please specify serial number				E04159							
Software conforming to FDA "Celsius FDA Edition". Meets the requirements for the use of electronically stored data sets and electronic signatures as laid down in Regulation 21 CFR Part 11 of the US Food and Drug Administration (FDA). Base licence for the control of one unit				E05019							
Integration per additional unit (up to max.15 units) into an already existent FDA-software licence (E05019)				FDAQ4							
IQ check list with device-specific works test data as support for validation by customer				D00103							
OQ check list with device-specific works test data for one free-selectable temperature value, incl. temperature distribution survey at Memmert for 27 measuring points to DIN 12880:2007-05 as support for validation by customer. Validation at customer site on demand	D00104										
OQ check list with device-specific works test data for one free-selectable temperature value, incl. temperature distribution survey at Memmert for 5 measuring points to DIN 12880:2007-05 as support for validation by customer valid for one thermoshelf; ditto for further thermoshelves VO on demand (VO and VOcool only). Validation at customer site on demand				D00117							
OQ check list with device-specific works test data for one free-selectable humidity and temperature value, incl. temperature distribution survey at Memmert for 27 measuring points to DIN 12880:2007-05 as support for validation by customer (models HCP and CTC). Validation at customer site on demand				D00104							
External measuring instrument with sensors for daylight and UV-light, with additional measuring head for temperature and humidity. Product information on demand (models HCP)				B04714							

MODEL VARIANTS

9.60	
SingleDISPLAY	TwinDISPLAY
ControlCOCKPIT with one TFT display	ControlCOCKPIT with two TFT displays
AVAILABLE APPLIANCES	AVAILABLE APPLIANCES
UN/UNm / UF/UFm / IN/INm / IF/IFm / SN / SF / IPP / IPS	UNplus/UNmplus / UFplus/UFmplus / UF TS / UNpa INplus/INmplus / IFplus/IFmplus / SNplus / SFplus ICOmed / IPPplus / ICP / HPP / ICH
One high-resolution TFT colour display with touch-sensitive buttons for selection of functions	Two high-resolution TFT colour displays with touch-sensitive buttons for selection of functions
Available parameters on the ControlCOCKPIT: Temperature (Celsius or Fahrenheit), fan speed, exhaust air flap position, programme time	Available parameters on the ControlCOCKPIT: Temperature (Celsius or Fahrenheit), fan speed, exhaust air flap position, programme time, relative humidity, illumination, ${\rm CO_2}$
One temperature sensor Pt100 DIN class A in a 4-wire circuit	Two Pt100 sensors DIN class A in a 4-wire circuit for mutual monitoring, taking over functions in case of an error
	HeatBALANCE function for application specific adjustment of heat output distribution (balance) between the upper and lower heating groups in an adjustment range between -50 % and +50 % (not valid for models 30, HPP110, IPP110plus, ICP, ICH)
AtmoCONTROL software for reading out, managing and organising the data logger via Ethernet interface (temporary trial version can be downloaded). USB stick with AtmoCONTROL software available as accessory (on demand)	AtmoCONTROL software on a USB stick for programming, managing and transferring programmes via Ethernet interface or USB port
	ControlCOCKPIT with USB port for uploading programmes, reading out protocol logs, activating the User-ID function
	Displaying of already logged protocol data on the ControlCOCKPIT (max 10,000 values correspond to approx. 1 week)
Ethernet interface on the rear of the appliance for reading out the protocol log and for online logging	Ethernet interface on the rear of the appliance for reading out the protocol log and for uploading programmes and for online logging
Double overtemperature protection: Electronic temperature monitoring with freely adjustable monitoring temperature, for models U, I, S with option A6 TWW/TWB (protection class 3.1 or 2), mechanical temperature limiter TB acc. to DIN 12880	Multiple overtemperature protection: Electronic temperature monitoring TWW/TWB (protection class 3.1 or 2 resp. 3.3 for units with active cooling) and mechanical temperature limiter TB (protection class 1) acc. to DIN 12880, AutoSAFETY automatically adjusts to the set value within a freely adjustable tolerance range. Setting individual MIN / MAX values for over/undertemperature alarm and also for all other parameters such as relative humidity, CO ₂
PID microprocessor control with	integrated auto-diagnostic system
Structured stainless steel housing scratch-resis	stant robust and durable: rear of zinc-plated steel

Structured stainless steel housing, scratch-resistant, robust and durable; rear of zinc-plated steel

High-temperature connectors on the rear of the appliance for single-phase power connection according to country specific systems and IEC standards

Internal data logger with a storage capacity of at least 10 years

German, English, French, Spanish, Polish, Czech, Hungarian language settings available on the ControlCOCKPIT

Digital backwards counter with target time setting, adjustable from 1 minute to 99 days

The SetpointWAIT function guarantees that the process time does not start until the set temperature is reached at all measuring points – optional for temperature values recorded by the freely positionable Pt100 sensors inside the chamber

Adjustment of three calibration values for temperature and additional appliance specific parameters directly at the ControlCOCKPIT





HEATING AND DRYING OVENS

UNIVERSAL OVEN U

PASS-THROUGH OVEN UF TS

PARAFFIN OVEN UNpa

STERILISER S

VACUUM OVEN VC

COOLED VACUUM OVEN VOcoo

INCUBATORS

INCUBATOR I

CO, INCUBATOR ICOmed

COMPRESSOR-COOLED INCUBATOR ICP

PELTIER-COOLED INCUBATOR IPP

COOLED STORAGE INCUBATOR IPS

CLIMATE CHAMBERS

CONSTANT CLIMATE CHAMBER HPE

HIIMIDITY CHAMBER HCP

LIMATE CHAMBER ICH

ENVIRONMENTAL TEST CHAMBER CTC/TTC

WATERBATHS / OILBATHS

WATERBATH W

OILBATH O

Memmert GmbH + Co. KG P.O. Box 1720 | D-91107 Schwabach Tel. +49 9122 925-0 | Fax +49 9122 14585

E-Mail: sales@memmert.com facebook.com/memmert.family

The platform for experts: www.atmosafe.net

YOUR MEMMERT PARTNER